




LABOMED® 
ideas for vision





PRIMA DNT Mu Microscope

Binocular Tube: 0-210° tiltable head

Objective: f=250 mm with objective lens protector. Optional: NuVar system with variable focal length adjustment from 300-400 mm with fine focus mechanism, WD 220-332mm, variable focal length from 287-487, WD 200-400mm

Magnichanger: 3 step (0.6X, 1.0X, 1.6X).

Magnification: 10X Eyepieces (4.1X, 6.8X, 10.9X) 12.5X Eyepieces (5.1X, 8.5X, 13.6X)

Vertical Movement of Arm: 500 mm, power saving feature with light cut-off at park position

Built-in Filters: Yellow, Green

IPD Adjustment: 50-75mm


Diopter Adjustment: ± 5mm with diopter lock

Light Source: 50 watt LED with 80K LUX

Optional Accessories: Beam Splitters (Straight, Inclined), Range of Camera Adaptors, 6MP HD Camera with 30 fps (@1080p), Rotoplate, Double Iris, Assistant Binoculars



Contact: 90380 81446




Double Beam Splitter Inclined



Camera Adaptors



Double Iris Diaphragm



Rotoplate